

CUSTOMISE THIS INSTRUMENT WITH YOUR BRAND COLOUR AND LOGO

25, 50ml Plastic Semi-Automatic Injector

MODEL	VERSION	DOSE RANGE	INCREMENTS	CODE
25ml	Polysulfane Barrel and Metal Pushrod	0.5 – 2.5ml	0.5ml	PAS802
25ml	Polypropylene Barrel and Metal Pushrod	0.5 – 2.5ml	0.5ml	PAS1356
25ml	Polypropylene Barrel and Plastic Pushrod	0.5 – 2.5ml	0.5ml	PAS1337
50ml	Polysulfane Cylinder and Metal Pushrod	1.0 – 5.0ml	1.0ml	PAS803
50ml	Polypropylene Barrel and Metal Pushrod	1.0 – 5.0ml	1.0ml	PAS1357
50ml	Polypropylene Barrel and Plastic Pushrod	1.0 – 5.0ml	1.0ml	PAS1338

- Comfortable handle.
- Easy choice of 5 fixed dose settings.
- Cost effective, reliable and durable to save you money.
- Built in needle magazine saves you time.
- Ergonomic fingertip action reduces hand strain.
- Service Kits available.



INTERCHANGEABLE COLOURED PUSH KNOBS
FOR EASY PRODUCT IDENTIFICATION
WX1421



BUILT IN NEEDLE
MAGAZINE

