

CUSTOMISE THIS INSTRUMENT WITH YOUR BRAND COLOUR AND LOGO

25ml, 50ml Powermaster

- Fast, quick refill action.
- Powered by Liquid Propane Gas or an Air Compressor.
- Strong, lightweight construction.
- Ergonomic; minimal hand strain.
- Multiple use with interchangeable nozzle assemblies.

PORTABLE MODEL	DOSE RANGE	INCREMENTS	CODE
25ml Fixed Hook	2.5ml - 25ml	2.5ml	P25P-FH
25ml Pour-On	2.5ml - 25ml	2.5ml	P25P-PO
25ml Sheep Oral	2.5ml - 25ml	2.5ml	P25P-SO
50ml Cattle Oral	5ml - 50ml	2.5ml	P50P-CO
50ml Fixed Hook	5ml - 50ml	2.5ml	P50P-FH
50ml Pour-On Bent Nozzle	5ml - 50ml	2.5ml	P50P-POB
50ml Sheep Oral	5ml - 50ml	2.5ml	P50P-SO

PORTABLE MODEL	DOSE RANGE	INCREMENTS	CODE
25ml Fixed Hook	2.5ml - 25ml	2.5ml	P25A-FH
25ml Pour-On	2.5ml - 25ml	2.5ml	P25A-PO
25ml Sheep Oral	2.5ml - 25ml	2.5ml	P25A-SO
50ml Cattle Oral	5ml - 50ml	2.5ml	P50A-CO
50ml Fixed Hook	5ml - 50ml	2.5ml	P50A-FH
50ml Pour-On Bent Nozzle	5ml - 50ml	2.5ml	P50A-POB
50ml Sheep Oral	5ml - 50ml	2.5ml	P50A-SO

PORTABLE MODEL	DOSE RANGE	INCREMENTS	CODE
25ml Fixed Hook	2.5ml - 25ml	2.5ml	P25G-FH
25ml Pour-On	2.5ml - 25ml	2.5ml	P25G-PO
25ml Sheep Oral	2.5ml - 25ml	2.5ml	P25G-SO
50ml Cattle Oral	5ml - 50ml	2.5ml	P50G-CO
50ml Fixed Hook	5ml - 50ml	2.5ml	P50G-FH
50ml Pour-On Bent Nozzle	5ml - 50ml	2.5ml	P50G-POB
50ml Sheep Oral	5ml - 50ml	2.5ml	P50G-SO

